

## **$J(^1\text{H}, ^1\text{H})$ -resolved spectrum**

\*\*\*\*Capitalized words are icons in the menu. The italic text you have to type\*\*\*\*

### **Get a good 1D NMR**

Set the cursor at  $\pm 1$ ppm to the outside peaks. Type *movesw*

*hom2dj*

set *ni=64*

*nt=8* or *16*

*sw1=xx*

Take a look at your 1D NMR. Wide peaks, vinyl, etc may have a width of more than 50Hz. This is the standard value in J-resolve (*sw1*). Be sure you increase it to 100 or 150Hz otherwise you get folded peaks.

*go*

Standard acquisition time 16-45 min depending on *ni*

!!! The splitting in 2<sup>nd</sup> dimension is due to 1H only. If you have F or any other atom that may split the signal you'll see in the first dimension. For instance, if you have one proton and one F splitting the signal you'll see a doublet in the 2<sup>nd</sup> dimension (1H splitting) and a doublet in 1<sup>st</sup> dimension ( 19F splitting) !!!

!!! The spectrum must be tilted!!!

### **Plot/save:**

1) Define your window:

*wc=xx*

*wc2=xx* (the *xx* value should be between 100 and 250 (Define the size of the window. If *wc=wc2* you'll get a square plot. I would recommend the value of *wc* around 100)

2) Plot: PROJ / HPROJ MAX / PLOT / CANCEL / MAIN /  
MENU / DISPLAY / PLOT / ALL CONTOURS / PAGE

After HPROJ MAX you'll get the projection. Adjust the height with the mouse.

In J Resolve you do not need both projections.

3) Save file as \*.eps

PROJ / HPROJ MAX / PLOT / VPROJ MAX / PLOT / CANCEL / MAIN  
MENU / DISPLAY / PLOT / ALL CONTOURS / *page('file name.eps')* /